


## ExpellPlus Filter Tips

### COMPATIBILITY CHART

For tips compatibility with other pipette brands, please contact [info@capp.dk](mailto:info@capp.dk)




Capp Bravo, ecopipette and CappAero single pipette	0.2-2µl 0.5-10µl	0.2-2µl 0.5-10µl	0.2-2µl 0.5-10µl	0.2-2µl 0.5-10µl	2-20µl	5-50µl 10-100µl
Capp electronic pipettes and CappAero multichannel	0.5-10µl 2-20µl	0.5-10µl 2-20µl	0.5-10µl 2-20µl	0.5-10µl 2-20µl	x	5-50µl 10-100µl
Biohit eLine, mLine, Proline Plus	0.1-3µl 0.5-10µl	0.1-3µl 0.5-10µl	0.1-3µl 0.5-10µl	0.1-3µl 0.5-10µl	2-20µl	2-20µl 10-100µl
Brand Transferpette, Transferpette S, Transferpette Electronic	0.5-10µl	0.5-10 µl	0.5-10µl	0.5-10µl	2-20µl	x
Eppendorf Reference, Research Plus	0.1-2.5µl 0.5-10µl	0.5-10µl	0.1-2.5µl 0.5-10µl	0.1-2.5µl 0.5-10µl	2-20µl	x
Finnpipette Digital	0.5-10µl	0.5-10 µl	0.5-10µl	0.5-10µl	2-20 µl	5-50µl
Gilson Pipetman	0.2-2µl 1-10µl	0.2-2µl 1-10µl	0.2-2µl 1-10µl	0.2-2µl 1-10µl	2-20µl	2-20µl 10-100µl
Rainin Classic, universal shaft (non LTS)	0.1-2µl 0.5-10µl	0.1-2µl 0.5-10µl	0.2-2µl 1-10µl	0.1-2µl 0.5-10µl	2-20 µl	2-20µl 10-100µl
VWR, Labnet, HTL mechanical pipettes	0.1-2µl 0.5-10µl	0.1-2µl 0.5-10µl	0.2-2µl 1-10µl	0.1-2µl 0.5-10µl	2-20 µl	2-20µl 10-100µl

# Expell & ExpellPlus Tips

## COMPATIBILITY CHART

For tips compatibility with other pipette brands, please contact [info@capp.dk](mailto:info@capp.dk)



5130070C 4130075C 5030070C 5030075C 5030080C 5030095C	5030090C 5131090C	5130100C 5130105C 5130110C	5130115C	5130130C 4130135C 5130135C 5130135C 5130140C 5130155C	5130150C 5131150C	5130120C 5130121C 5130122C	5130123C	5130160C 5130165C	5130190C 5130195C
2-20µl 5-50µl 10-100µl 20-200µl	2-20µl 5-50µl 10-100µl 20-200µl	2-20µl 5-50µl 10-100µl 20-200µl	2-20µl 5-50µl 10-100µl 20-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	1-5ml	1-10ml
5-50µl 10-100µl 20-200µl	5-50µl 10-100µl 20-200µl	5-50µl 10-100µl 20-200µl 30-300µl	5-50µl 10-100µl 20-200µl 30-300µl	100-1200µl	100-1200µl	100-1200µl	100-1200µl	x	x
2-20µl 10-100µl 20-200µl	10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	x	1-10ml
2-20µl 10-100µl 20-200µl	20-200µl	20-200µl 15-300µl	20-200µl 15-300µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	x	x
2-20µl 10-100µl 20-200µl	20-200µl	10-100µl 20-200µl 30-300µl	10-100µl 20-200µl 30-300µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	x	1-10ml
2-20µl 5-50µl 20-200µl	5-50µl 20-200µl	5-50µl 20-200µl 30-300µl	5-50µl 20-200µl 30-300µl	100-1000µl 200-1000µl	100-1000µl 200-1000µl	100-1000µl 200-1000µl	100-1000µl 200-1000µl	x	x
2-20µl 20-100µl 50-200µl	2-20µl 10-100µl 20-200µl	2-20µl 20-100µl 50-200µl	2-20µl 20-100µl 50-200µl	100-1000µl 200-1000µl	100-1000µl 200-1000µl	100-1000µl 200-1000µl	100-1000µl 200-1000µl	1-5ml	1-10ml
2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	1-5ml	1-10ml
2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	2-20µl 10-100µl 20-200µl	100-1000µl	100-1000µl	100-1000µl	100-1000µl	1-5ml	1-10ml